

Product Datasheet

Human Cartpt ELISA Kit (orb406696)

Description	Human Cocaine-and amphetamine-regulated transcript protein (Cartpt) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	Cartpt
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q16568
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.