

## Product Datasheet

### Human CD28 ELISA Kit (orb406702)

|                          |                                                             |
|--------------------------|-------------------------------------------------------------|
| <b>Description</b>       | human soluble cluster of differentiation 28,sCD28 ELISA Kit |
| <b>Reactivity</b>        | Human                                                       |
| <b>Range</b>             | 0.625 ng/mL-40 ng/mL                                        |
| <b>Target</b>            | CD28                                                        |
| <b>Note</b>              | For research use only                                       |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                           |
| <b>Uniprot ID</b>        | <b>P10747</b>                                               |
| <b>Sensitivity</b>       | 0.156 ng/mL                                                 |
| <b>Expiration Date</b>   | Please enquire.                                             |