

Product Datasheet

Human sAPP α ELISA Kit (orb406714)

Description	Human soluble amyloid precursor protein α (sAPP α) ELISA kit
Reactivity	Human
Range	12.5 ng/mL-800 ng/mL
Target	sAPP α
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates, cerebrospinal fluid (CSF)
Sensitivity	3.12 ng/mL
Expiration Date	Please enquire.