

Product Datasheet

Human RANKL ELISA Kit (orb406717)

Description	Human soluble receptor activator of nuclear factor-kB ligand,sRANKL ELISA Kit
Reactivity	Human
Range	7.8 pg/mL-500 pg/mL
Target	TNFSF11
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supeernates, tissue homogenates
Uniprot ID	O14788
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.