

Product Datasheet

Human sLOX-1 ELISA Kit (orb406720)

Description	Human soluble lectin-like oxidized low density lipoprotein receptor-1(sLOX-1)ELISA Kit
Reactivity	Human
Range	2.74 pg/ ml - 2000 pg /ml
Target	sLOX-1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates.
Uniprot ID	P78380
Sensitivity	2.74 pg/ ml
Expiration Date	Please enquire.