

Product Datasheet

Human sPECAM1/sCD31 ELISA Kit (orb406724)

Description	Human soluble platelet endothelial cell adhesion molecule 1,sPECAM-1 ELISA Kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	sPECAM1/sCD31
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates
Uniprot ID	P16284
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.