

Product Datasheet

Human TNFSF10 ELISA Kit (orb406726)

Description	Human soluble tumor necrosis factor-related apoptosis inducing ligand,sTRAIL ELISA KIT
Reactivity	Human
Range	15.6 pg/ml - 1000pg/ml
Target	TNFSF10
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P50591
Sensitivity	3.9 pg/ml
Expiration Date	Please enquire.