

Product Datasheet

Human NOTCH3 ELISA Kit (orb406734)

Description	Human Neurogenic locus notch homolog protein 3(NOTCH3) ELISA kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Target	NOTCH3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9UM47
Sensitivity	31.25 pg/mL
Expiration Date	Please enquire.