

Product Datasheet

Human NQO1 ELISA Kit (orb406735)

Description	Human NAD (P)H Dehydrogenase, Quinone 1 (NQO1)ELISA Kit
Reactivity	Human
Range	39 pg/mL-2500 pg/mL
Target	NQO1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P15559
Sensitivity	9.75 pg/mL
Expiration Date	Please enquire.