

Product Datasheet

Human KDM5A ELISA Kit (orb406736)

Description	Human Lysine-specific demethylase 5A(KDM5A) ELISA kit
Reactivity	Human
Range	23.5 pg/mL-1500 pg/mL
Target	KDM5A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P29375
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.