

Product Datasheet

Human TDP1 ELISA Kit (orb406742)

Description	Human Tyrosyl-DNA Phosphodiesterase 1(TDP1)ELISA Kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	TDP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9NUW8
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.