

Product Datasheet

Human HLA-B27 ELISA Kit (orb406744)

Description	Human leucocyte antigen B27,HLA-B27 ELISA Kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	HLA-B27
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, urine, cerebrospinal fluid (CSF)
Uniprot ID	P03989
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.