

Product Datasheet

Human HLA-E ELISA Kit (orb406745)

Description	Human major histocompatibility complex, class I, E (HLA-E) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	HLA-E
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P13747
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.