

Product Datasheet

Human pMAPT /pTAU ELISA Kit (orb406767)

Description	Human phosphorylated microtubule-associated protein tau (pMAPT /pTAU) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	pMAPT /pTAU
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cerebrospinal fluid (CSF)
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.