

Product Datasheet

Human GPLD1 ELISA Kit (orb406779)

Description	Human Phosphatidylinositol-glycan-specific phospholipase D(GPLD1) ELISA kit
Reactivity	Human
Range	3.12 ng/mL-200 ng/mL
Target	GPLD1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P80108
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.