

Product Datasheet

Human TPK1 ELISA Kit (orb406794)

Description	Human Thiamin pyrophosphokinase 1 (TPK1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	TPK1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9H3S4
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.