

Product Datasheet

Human SOD2 ELISA Kit (orb406823)

Description	Human Superoxide dismutase [Mn], mitochondrial (SOD2) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	SOD2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Uniprot ID	P04179
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.