

## Product Datasheet

### Human ANXA2 ELISA Kit (orb406838)

|                          |                                       |
|--------------------------|---------------------------------------|
| <b>Description</b>       | Human Annexin II(ANX-II)<br>ELISA Kit |
| <b>Reactivity</b>        | Human                                 |
| <b>Range</b>             | 0.39 ng/mL-25 ng/mL                   |
| <b>Target</b>            | ANXA2                                 |
| <b>Note</b>              | For research use only                 |
| <b>Application notes</b> | Detection Wavelength: 450 nm          |
| <b>Sample Types</b>      | serum, urine, saliva, cell<br>lysates |
| <b>Uniprot ID</b>        | <b>P07355</b>                         |
| <b>Sensitivity</b>       | 0.098 ng/mL                           |
| <b>Expiration Date</b>   | Please enquire.                       |