

Product Datasheet

Human Annexin A3 ELISA Kit (orb406839)

Description	Human Annexin III(ANX-III) ELISA Kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Target	ANXA3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P12429
Sensitivity	0.195 ng/mL
Expiration Date	Please enquire.