

Product Datasheet

Human Autotaxin ELISA Kit (orb406847)

Description	Human Ectonucleotide pyrophosphatase/phosphodiesterase family member 2 (ENPP2) ELISA kit
Reactivity	Human
Range	25 ng/mL-2000 ng/mL
Target	ENPP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q13822
Sensitivity	31.25 ng/mL
Expiration Date	Please enquire.