

Product Datasheet

Human ENPP3 ELISA Kit (orb406848)

Description	Human Ectonucleotide pyrophosphatase/phosphodiesterase family member 3(ENPP3) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	ENPP3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O14638
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.