

Product Datasheet

Human PGRMC1 ELISA Kit (orb406849)

Description	Human Membrane-associated progesterone receptor component 1(PGRMC1) ELISA kit
Reactivity	Human
Range	47 pg/mL-3000 pg/mL
Target	PGRMC1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O00264
Sensitivity	11.7 pg/mL
Expiration Date	Please enquire.