

Product Datasheet

Human ESM1 ELISA Kit (orb406861)

Description	Human endothelial cell-specific molecule 1(ESM-1) ELISA Kit
Reactivity	Human
Range	3.12 ng/ml-200 ng/ml.
Target	ESM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenatess.
Uniprot ID	Q9NQ30
Sensitivity	0.78 ng/ml
Expiration Date	Please enquire.