

Product Datasheet

Human OTC ELISA Kit (orb406870)

Description	Human Ornithine carbamoyltransferase, mitochondrial(OTC) ELISA kit
Reactivity	Human
Range	1.25 mU/mL-80 mU/mL
Target	OTC
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P00480
Sensitivity	0.312 mU/mL
Expiration Date	Please enquire.