

Product Datasheet

Human PLAUR ELISA Kit (orb406875)

Description	Human Urokinase-type Plasminogen Activator Receptor, PLAUR/uPAR ELISA kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	PLAUR
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q03405
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.