

Product Datasheet

Human PIVKA-II ELISA Kit (orb406883)

Description	Human Protein Induced Vitamin K Absence or Antagonist-II(PIVKA-II) ELISA Kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	PIVKA-II
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, ascitic fluid, tissue homogenates
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.