

## Product Datasheet

### Human CORT ELISA Kit (orb406900)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Human Corticosterone,CORT<br>ELISA Kit |
| <b>Reactivity</b>        | Human                                  |
| <b>Range</b>             | 0.5 ng/mL-20 ng/mL                     |
| <b>Target</b>            | CORT                                   |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Detection Wavelength: 450 nm           |
| <b>Sample Types</b>      | serum, plasma                          |
| <b>Sensitivity</b>       | 0.5 ng/mL                              |
| <b>Expiration Date</b>   | Please enquire.                        |