

Product Datasheet

Human PIM1 ELISA Kit (orb406914)

Description	Human Proto-oncogene serine/threonine-protein kinase pim-1(PIM1) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	PIM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	P11309
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.