

## Product Datasheet

### Human dynorphin ELISA Kit (orb406934)

|                          |                                   |
|--------------------------|-----------------------------------|
| <b>Description</b>       | Human dynorphin,Dyn ELISA Kit     |
| <b>Reactivity</b>        | Human                             |
| <b>Range</b>             | 3.12 pg/mL-200 pg/mL              |
| <b>Target</b>            | Dyn                               |
| <b>Note</b>              | For research use only             |
| <b>Application notes</b> | Detection Wavelength: 450 nm      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates |
| <b>Sensitivity</b>       | 0.78 pg/mL                        |
| <b>Expiration Date</b>   | Please enquire.                   |