

Product Datasheet

Human SIRT1 ELISA Kit (orb406957)

Description	Human NAD-dependent deacetylase sirtuin-1 (SIRT1/SIR2L1) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	SIRT1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q96EB6
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.