

Product Datasheet

CAMK1D antibody (orb40696)

Description	CAMK1D antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Calcium/calmodulin-dependent protein kinase type 1D protein (1-180AA)
Target	CAMK1D
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IU85
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.