

## Product Datasheet

### Human HSP-70 ELISA Kit (orb406972)

|                          |                                                     |
|--------------------------|-----------------------------------------------------|
| <b>Description</b>       | Human heat shock protein 70(HSP-70) (IgG) ELISA kit |
| <b>Reactivity</b>        | Human                                               |
| <b>Range</b>             | Request Information                                 |
| <b>Target</b>            | HSP70 Ab                                            |
| <b>Note</b>              | For research use only                               |
| <b>Application notes</b> | Detection Wavelength: 450 nm                        |
| <b>Sample Types</b>      | serum, plasma                                       |
| <b>Sensitivity</b>       | Request Information                                 |
| <b>Expiration Date</b>   | Please enquire.                                     |