

Product Datasheet

Human HSP-72 ELISA Kit (orb406973)

Description	Human Heat Shock Protein 72(HSP-72) ELISA Kit
Reactivity	Human
Range	3.12 ng/mL-200 ng/mL
Target	HSP72
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, cell lysates
Uniprot ID	P0DMV8
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.