

Product Datasheet

Human TRAP1 ELISA Kit (orb406974)

Description	Human Heat shock protein 75 kDa, mitochondrial (TRAP1/HSP75) ELISA kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	TRAP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q12931
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.