

Product Datasheet

Human PSG3 ELISA Kit (orb406982)

Description	Human Pregnancy-specific beta-1-glycoprotein 3(PSG3) ELISA kit
Reactivity	Human
Range	0.625 ng/mL-50 ng/mL
Target	PSG3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q16557
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.