

Product Datasheet

Human HLAMP-2 ELISA Kit (orb406990)

Description	Human lysosomal-associated membrane protein 2,HLAMP-2 ELISA Kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	LAMP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P13473
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.