

Product Datasheet

Human NCAM1 ELISA Kit (orb407051)

Description	Human Neural cell adhesion molecule 1 (NCAM1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	NCAM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P13591
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.