

Product Datasheet

Human C10orf59 ELISA Kit (orb407059)

Description Human Renalase/ Monoamine Oxidase-C(RNLS/ MAO-C)ELISA Kit

Reactivity Human

Range 0.625 µg/mL-40 µg/mL

Target C10orf59

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates

Uniprot ID **Q5VYX0**

Sensitivity 0.156 µg/mL

Expiration Date Please enquire.