

## Product Datasheet

### Human OX-B ELISA Kit (orb407078)

|                          |                                      |
|--------------------------|--------------------------------------|
| <b>Description</b>       | Human Orexin B,OX-B ELISA Kit        |
| <b>Reactivity</b>        | Human                                |
| <b>Range</b>             | 62.5 pg/mL-4000 pg/mL                |
| <b>Target</b>            | OX-B                                 |
| <b>Note</b>              | For research use only                |
| <b>Application notes</b> | Detection Wavelength: 450 nm         |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates |
| <b>Sensitivity</b>       | 15.6 pg/mL                           |
| <b>Expiration Date</b>   | Please enquire.                      |