

Product Datasheet

Human TRPV1 ELISA Kit (orb407099)

Description	Human transient receptor potential cation channel subfamily V member 1 (TrpV1)ELISA Kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	TRPV1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8NER1
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.