

Product Datasheet

Human LMTK3 ELISA Kit (orb407102)

Description	Human Serine/threonine-protein kinase LMTK3(LMTK3) ELISA kit
Reactivity	Human
Range	18.75 pg/mL-1200 pg/mL
Target	LMTK3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q96Q04
Sensitivity	4.7 pg/mL
Expiration Date	Please enquire.