

Product Datasheet

Human SGK1 ELISA Kit (orb407104)

Description	Human Serine/threonine-protein kinase Sgk1(SGK1) ELISA kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Target	SGK1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O00141
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.