

Product Datasheet

Human SPINK6 ELISA Kit (orb407108)

Description	Human Serine protease inhibitor Kazal-type 6(SPINK6) ELISA kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	SPINK6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q6UWN8
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.