

Product Datasheet

Human MASP2 ELISA Kit (orb407109)

Description	Human mannan-binding lectin serine peptidase 2 (MASP2)ELISA kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Target	MASP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, cerebrospinal fluid (CSF)
Uniprot ID	O00187
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.