

## Product Datasheet

### Human FGF1 ELISA Kit (orb407119)

|                          |                                                            |
|--------------------------|------------------------------------------------------------|
| <b>Description</b>       | Human acidic fibroblast growth factor,aFGF/FGF-1 ELISA Kit |
| <b>Reactivity</b>        | Human                                                      |
| <b>Range</b>             | 15.6 pg/mL-1000 pg/mL                                      |
| <b>Target</b>            | FGF1                                                       |
| <b>Note</b>              | For research use only                                      |
| <b>Application notes</b> | Detection Wavelength: 450 nm                               |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                          |
| <b>Uniprot ID</b>        | <b>P05230</b>                                              |
| <b>Sensitivity</b>       | 3.9 pg/mL                                                  |
| <b>Expiration Date</b>   | Please enquire.                                            |