

Product Datasheet

Human SMPD1 ELISA Kit (orb407120)

Description	Human acid sphingomyelinase, ASM ELISA Kit
Reactivity	Human
Range	1.25 ng/mL-80 ng/mL
Target	SMPD1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, cerebrospinal fluid (CSF), urine
Uniprot ID	P17405
Sensitivity	0.31 ng/mL
Expiration Date	Please enquire.