

## Product Datasheet

### Human MBP ELISA Kit (orb407122)

|                          |                                             |
|--------------------------|---------------------------------------------|
| <b>Description</b>       | Human myelin basic protein,MBP<br>ELISA Kit |
| <b>Reactivity</b>        | Human                                       |
| <b>Range</b>             | 7.8 ng/ml- 500ng/ml                         |
| <b>Target</b>            | MBP                                         |
| <b>Note</b>              | For research use only                       |
| <b>Application notes</b> | Detection Wavelength: 450 nm                |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates        |
| <b>Uniprot ID</b>        | <b>P02686</b>                               |
| <b>Sensitivity</b>       | 1.95 ng/ml                                  |
| <b>Expiration Date</b>   | Please enquire.                             |