

Product Datasheet

Human TREM1 ELISA Kit (orb407127)

Description	Human Triggering Receptor Expresses on Myeloid Cells-1,TREM-1 ELISA Kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	TREM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9NP99
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.