



Product Datasheet

CSNK1G2 antibody (orb40714)


Description	CSNK1G2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human Casein kinase I isoform gamma-2 protein (316-415AA)
Target	CSNK1G2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P78368
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human testis tissue.



Immunohistochemical staining of human brain tissue.



Western blot analysis of Mouse gonadal tissue.