

Product Datasheet

Human Ferritin ELISA Kit (orb407170)

Description	Human Ferritin, Heavy Polypeptide (FTH) ELISA Kit
Reactivity	Human
Range	0.45 ng/mL-30 ng/mL
Target	FTH1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P02794
Sensitivity	0.11 ng/mL
Expiration Date	Please enquire.